

Deltalab Large Nitrile Gloves (Powder Free)

Product Code: 402NG

Brand: Deltalab

Suitable for medical and laboratory use, Deltalab's Large Nitrile Gloves (Powder Free) are made of blue nitrile.

These maximum safety and waterproof gloves are three times more than resistant than latex gloves. With a microtextured outer surface on the fingertips, they provide greater grip strength. Their chlorinated inner surface is smooth and comfortable.

Made of materials suitable for food use and free of protein, they can be used with a temperature range of -20°C to 50°C and can be easily and conveniently removed.



Specifications:

Thickness	Finger	0.10mm (± 0.02mm)
	Palm	0.07mm (± 0.02mm)
	Fist	0.05mm (± 0.02mm)
Expiration	60 months from manufacturing date	
Size	Large	
Material	Nitrile	
Sterility	Non-Sterile	
Quantity	100	

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com