

Deltalab Microtitre Plates (F Well)

Product Code: 90001 | SP

Brand: Deltalab

Deltalab's Microtitre Plates (F Well) are standard 96 well-plates made of crystal-clear polystyrene.

They are suitable for various applications such as diluting, sample preparation, and protein precipitation. These flat bottom plates can also be used for cell growth in suspension and for ELISA applications.

Features:

- Fully compatible to all current standard-equipment (reading and washing devices).
- Can be safely stacked, causing no danger of spillage.
- Appropriate gripping surfaces at the side.
- Frosted surface of the plates prevent errors from extraneous light.
- Alpha-numerical for easy sample location.
- Space for identification marking on the side surface.
- Covers with condensation-rings to avoid medium-losses and uneven evaporation.
- Lid orientation corners reduce risk of cross-contamination.

Specifications:

Plate Dimensions	127.7 x 85.8mm
Well Plate Estimate Value	0.404µl
Material	Polystyrene
Sterility	Non-Sterile
Quantity	100

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

