


## Detalab Non-Sterile General Purpose Pipettes (8ml)

**Product Code:** 200000

**Brand:** Detalab

Made of transparent polyethylene, Detalab's Non-Sterile General Purpose Pipettes (8ml) can easily transfer liquids without the risk of cross-contamination occurring. Their cohesive design with a built-in pipette bulb makes them safe and comfortable.



### Specifications:

<b>Volumetrical Tolerance</b>	± 10%
<b>Material</b>	Polyethylene
<b>Sterility</b>	Non-Sterile
<b>Length (mm)</b>	150
<b>Stem (Ø mm)</b>	6.3
<b>Total Capacity (ml)</b>	8
<b>Bulb Draw (ml)</b>	4.6
<b>Quantity</b>	400

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)