

## Deltalab Non-Sterile Graduated Pipettes (1ml)

**Product Code:** 200003

**Brand:** Deltalab

Deltalab's Non-Sterile Graduated Pipettes (1ml) are made of transparent, inert, and low density polyethylene. Suitable for laboratories in various sectors, these pasteur pipettes can be used for sampling or measuring transfers of aqueous solutions.



### Specifications:

<b>Volumetrical Tolerance</b>	±10%
<b>Presentation</b>	6 x 500
<b>Sterility</b>	Non-Sterile
<b>Length (mm)</b>	155
<b>Stem (Ø mm)</b>	4.6
<b>Graduated Until (ml)</b>	1
<b>Total Capacity (ml)</b>	5.8
<b>Bulb Draw (ml)</b>	3.4
<b>Graduations (ml)</b>	0.25
<b>Material</b>	Polyethylene

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)