


Deltalab Non-Sterile Graduated Pipettes (3ml)

Product Code: 200006

Brand: Deltalab

Deltalab's Non-Sterile Graduated Pipettes (3ml) are made of transparent, inert, and low density polyethylene. Suitable for laboratories in various sectors, these pasteur pipettes can be used for sampling or measuring transfers of aqueous solutions.



Specifications:

Volumetrical Tolerance	±10%
Presentation	10 x 500
Sterility	Non-Sterile
Length (mm)	155
Stem (Ø mm)	7.8
Graduated Until (ml)	3
Total Capacity (ml)	7.5
Bulb Draw (ml)	3.2
Graduations (ml)	0.50
Material	Polyethylene

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com