


## Deltalab Non-Sterile Micro Tip Pasteur Pipettes (1ml)

**Product Code:** 210003

**Brand:** Deltalab

Made of transparent and inert polyethylene, Deltalab's Non-Sterile Micro Tip Pasteur Pipettes (1ml) can be used for fragile, small drop work.



### Specifications:

<b>Volumetrical Tolerance</b>	± 10%
<b>Material</b>	Polyethylene
<b>Sterility</b>	Non-Sterile
<b>Length (mm)</b>	104
<b>Stem (Ø mm)</b>	3
<b>Total Capacity (ml)</b>	1.5
<b>Bulb Draw (ml)</b>	1
<b>Quantity</b>	400

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)