

## Deltalab Non-Sterile Microtitre Plates (LI Well)

**Product Code:** 900011.10

**Brand:** Deltalab

Deltalab's Non-Sterile Microtitre Plates (LI Well) are standard 96 well-plates made of crystal-clear polystyrene.


They are suitable for various applications such as diluting, sample preparation, and protein precipitation. These flat bottom plates can also be used for cell growth in suspension and for ELISA applications.

### Features

- Fully compatible to all current standard-equipment (reading and washing devices)
- Can be safely stacked, causing no danger of spillage
- Appropriate gripping surfaces at the side
- Frosted surface of the plates prevent errors from extraneous light
- Alpha-numerical for easy sample location
- Space for identification marking on the side surface
- Covers with condensation-rings to avoid medium-losses and uneven evaporation
- Lid orientation corners reduce risk of cross-contamination

### Specifications

<b>Plate Dimensions (mm)</b>	127.7 x 85.8
<b>Material</b>	Polystyrene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	50



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)