

Deltalab Polyethylene Wide Mouth Jars with Lid & Insert Plug (250ml)

Product Code: 202802

Brand: Deltalab

The high density Deltalab Polyethylene Wide Mouth Jars with Lid & Insert Plug (250ml) are useful for storing dry chemicals.

The jars and plugs are naturally translucent and the lids are black.

Upon request, the jars can be supplied individually wrapped for a minimum order of 3 cases for each code, and sterile for a minimum order of 6 cases, for each code.



Specifications:

| | |
|-----------------------------|--------------|
| Capacity (ml) | 250 |
| Internal Mouth (mm) | 54 |
| Body (mm) | 67 |
| Height with Lid (mm) | 100 |
| Alimentary Use? | Yes |
| Material | Polyethylene |
| Sterility | Non-Sterile |
| Quantity | 125 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com