

## Deltalab Polyethylene Wide Mouth Jars with Lid & Insert Plug (500ml)

**Product Code:** 202803

**Brand:** Deltalab

The high density Deltalab Polyethylene Wide Mouth Jars with Lid & Insert Plug (500ml) are useful for storing dry chemicals.

The jars and plugs are naturally translucent and the lids are black.

Upon request, the jars can be supplied individually wrapped for a minimum order of 3 cases for each code, and sterile for a minimum order of 6 cases, for each code.



### Specifications:

<b>Capacity (ml)</b>	500
<b>Internal Mouth (mm)</b>	67
<b>Body (mm)</b>	80
<b>Height with Lid (mm)</b>	131
<b>Alimentary Use?</b>	Yes
<b>Material</b>	Polyethylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	125

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)