


Deltalab Radiation Sterilised Graduated Pipettes (3ml)

Product Code: 200037

Brand: Deltalab

Made of transparent and inert low density polyethylene, Deltalab's Radiation Sterilised Graduated Pipettes (3ml) are ideal for general laboratory procedures that need samples to be transferred safely and quickly.



Specifications:

Volumetrical Tolerance	± 10%
Presentation	6 x 500 individually-wrapped flow pack
Sterility	Radiation Sterilised
Length (mm)	154
Stem (Ø mm)	7.8
Graduated Until (ml)	3
Total Capacity (ml)	7
Bulb Draw (ml)	3.2
Graduations (ml)	0.50
Material	Polyethylene

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com