

Deltalab Rectangular Waste Containers (30 Litres)

Product Code: 240035

Brand: Deltalab

Made of perforated and resistant virgin polypropylene, Deltalab's Rectangular Waste Containers (30 Litres) are suitable for solid or semi-solid waste (II and III groups) and are manufactured following new French standard NFX 30 505.

Specifications:

Capacity (litres)	30
Container Dimensions (mm)	415 x 314 x 373
Maximum Load (kg)	14
Weight (ut. g)	1.4
Minimum Wall Thickness (mm)	2.5
Material	Polypropylene
Sterility	Non-Sterile
Quantity	1



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com