

Deltalab Round Safety Containers (25 Litres)

Product Code: 240028


Brand: Deltalab

Suitable for solid waste types II and III, Deltalab's Round Safety Containers (25 Litres) are manufactured in virgin raw materials enabling incineration and in accordance with ONU standards and major European standards.

They are compound by a body bearing the hazard warning label, a handle, and two lids; allowing the container to be totally opened, taking advantage of the total diameter of the mouth of the container and making it very useful to throw away bigger volume waste.

Specifications:

Capacity (litres)	25
Material	Polyethylene
Mouth Diameter (mm)	323
Base Diameter (mm)	283
Height (mm)	387
Sterility	Non-Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com