

Deltalab Round Safety Containers (7 Litres)

Product Code: 240007

Brand: Deltalab

Made of polyethylene and suitable for solid waste types II and III, Deltalab's Round Safety Containers (7 Litres) are manufactured with raw materials to enable incineration and are in accordance with major European and international standards (ONU ADR, etc.).

They are compound by a body bearing the hazard warning label, two lids, and a handle. The first lid offers an entire opening, taking advantage of the total diameter of the mouth of the container; very useful to throw away bigger volume waste. It can be closed hermetically and definitively. The second lid presents a partial opening, which may be closed temporarily or definitively. It embodies two devices in order to easy the extraction/disconnection on dental and analytical needles, and a 70 x 42mm rectangular orifice.

Specifications:

| | |
|----------------------------------|--------------|
| Capacity (litres) | 7 |
| Material (body & lid) | Polyethylene |
| Mouth Diameter (mm) | 210 |
| Base Diameter (mm) | 178 |
| Height (mm) | 240 |
| Sterility | Non-Sterile |
| Quantity | 30 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com