

Deltalab Security Screw Cap Containers (50ml)

Product Code: 202841

Brand: Deltalab

Deltalab's Security Screw Cap Containers (50ml) are ultra clear polypropylene containers with yellow polyethylene caps and moulded graduations in ml.

Supplied uncapped with the caps and containers in the same box, these laboratory containers are ideal for the pneumatic transport of liquid samples, including those that have phormaline, because they have a watertight closure. The cap, with an innovative shape, has a double closure. The base and the cap are ribbed so the handling with gloves is comfortable and reliable.

There is a ring in the upper part of the body in order to avoid the dripping in case of liquid decanting.



Specifications:

| | |
|--------------------------|---|
| Grad. Volume (ml) | 40 |
| Total Volume (ml) | 50 |
| Mouth (Ø mm) | 45 |
| Height (mm) | 43 |
| Material(s) | Polypropylene (containers) Polyethylene (caps) |
| Sterility | Non-Sterile |
| Quantity | 600 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com