

## Deltalab Security Screw Cap Containers (76ml)

**Product Code:** 202842

**Brand:** Deltalab

Deltalab's Security Screw Cap Containers (76ml) are ultra clear polypropylene containers with yellow polyethylene caps and moulded graduations in ml.

Supplied uncapped with the caps and containers in the same box, these laboratory containers are ideal for the pneumatic transport of liquid samples, including those that have phormaline, because they have a watertight closure. The cap, with an innovative shape, has a double closure. The base and the cap are ribbed so the handling with gloves is comfortable and reliable.

There is a ring in the upper part of the body in order to avoid the dripping in case of liquid decanting.



### Specifications:

<b>Grad. Volume (ml)</b>	60
<b>Total Volume (ml)</b>	76
<b>Mouth (Ø mm)</b>	45
<b>Height (mm)</b>	60
<b>Material(s)</b>	Polypropylene (containers) Polyethylene (caps)
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)