

## Deltalab Small Nitrile Gloves (Powder Free)

**Product Code:** 402NP

**Brand:** Deltalab

Suitable for medical and laboratory use, Deltalab's Small Nitrile Gloves (Powder Free) are made of blue nitrile.

These maximum safety and waterproof gloves are three times more resistant than latex gloves. With a microtextured outer surface on the fingertips, they provide greater grip strength. Their chlorinated inner surface is smooth and comfortable.

Made of materials suitable for food use and free of protein, they can be used with a temperature range of -20°C to 50°C and can be easily and conveniently removed.



### Specifications:

<b>Thickness</b>	<b>Finger</b>	0.10mm (± 0.02mm)
	<b>Palm</b>	0.07mm (± 0.02mm)
	<b>Fist</b>	0.05mm (± 0.02mm)
<b>Expiration</b>	60 months from manufacturing date	
<b>Size</b>	Small	
<b>Material</b>	Nitrile	
<b>Sterility</b>	Non-Sterile	
<b>Quantity</b>	100	

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)