

## Deltalab Snap Top Bucket (1000ml)

**Product Code:** 241013

**Brand:** Deltalab

Deltalab's Snap Top Buckets (1000ml) feature a body and self sealing security lid made of white autoclavable polypropylene.

Made with materials suitable for alimentary use, these laboratory containers are supplied uncapped.



### Specifications:

<b>Nominal Capacity (ml)</b>	1000
<b>Dimensions (Internal Mouth x Base x Height)</b>	107 x 94 x 144mm
<b>Weight (g)</b>	36
<b>Material</b>	Polypropylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)