

## Deltalab Square Petri Dishes (120 x 120mm)

**Product Code:** 200204

**Brand:** Deltalab

Made of polystyrene, Deltalab's Square Petri Dishes (120 x 120mm) come with four vents and are aseptic. These petri dishes are suitable for automatic filling machines.



### Specifications:

<b>Dimensions (mm)</b>	120 x 120
<b>Material</b>	Polystyrene
<b>Sterility</b>	Sterile
<b>Quantity</b>	240

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)