

Deltalab Sterile Amies Charcoal Swabs (Aluminium & Viscose)

Product Code: 300281/1

Brand: Deltalab

Deltalab's Sterile Amies Charcoal Swabs (Aluminium & Viscose) are a modification of Cary Blair media, which is at the same time a modification of Stuart media. It assures the viability of organisms such as *Trichomonas Vaginalis*, *Streptococcus Pyogenes*, and *Streptococcus Pneumoniae*.

Swabs are sterilised by radiation.
Peel-packed.

Specifications

| | |
|----------------------------------|---------------------------------------|
| Peel-Pack Dimensions (mm) | 38 x 210 |
| Material(s) | Aluminium and Viscose (with charcoal) |
| Expiry Date | 30 months after sterilisation date |
| Quantity | 100 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com