


Deltalab Sterile Micro Tip Pasteur Pipettes (1.5ml)

Product Code: 210023

Brand: Deltalab

Made of transparent and inert polyethylene, Deltalab's Sterile Micro Tip Pasteur Pipettes (1.5ml) can be used for fragile, small drop work.

Each pipette is sterilised by ethylene oxide.



Specifications:

Volumetrical Tolerance	± 10%
Material	Polyethylene
Sterility	Sterile
Length (mm)	104
Stem (Ø mm)	3
Total Capacity (ml)	1.5
Bulb Draw (ml)	1
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com