

Deltalab Sterile Water Sampling Jars - With Sodium Thiosulfate (500ml)

Product Code: 282340


Brand: Deltalab

Deltalab's Sterile Water Sampling Jars - With Sodium Thiosulfate (500ml) can be used to analyze water meant for human consumption, swimming pools, and any other water where chlorine may modify the composition of the sample while being transported.

Each bottle is labelled specifying description, code, lot, and expiry date.

Features:

- Transparent. Made of PET (Polyethylene Terephthalate).
- Wide mouth and grips at two sides, easing the sample collection.
- Red cap made of polyethylene with an internal liner made of polexan with a label that performs as tamper evident.



Specifications:

Internal Mouth Diameter (mm)	55
Capacity (ml)	500
Jar Weight (g)	44
Sterility	Sterilised by radiation
Material(s)	Polyethylene Terephthalate (jar) and Polyethylene (cap)
Quantity	44

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com