

## Deltalab Tamper Evident Bucket (2500ml)

**Product Code:** 222802

**Brand:** Deltalab

Deltalab's Tamper Evident Bucket (2500ml) is a disposable container made of autoclavable polypropylene. It features a tight fitting and leakproof lid.

It also has a sturdy and convenient white plastic handle for easier carrying.

It is supplied with the lid apart and made with materials suitable for alimentary use.



### Specifications:

|                                     |               |
|-------------------------------------|---------------|
| <b>Capacity (ml)</b>                | 2500          |
| <b>External Lid Diameter (mm)</b>   | 200           |
| <b>Internal Mouth Diameter (mm)</b> | 184           |
| <b>Height Without Lid (mm)</b>      | 138           |
| <b>Body Weight (g)</b>              | 105           |
| <b>Material</b>                     | Polypropylene |
| <b>Sterility</b>                    | Non-Sterile   |
| <b>Quantity</b>                     | 1             |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)