

Deltalab Tamper Evident Bucket (5000ml)

Product Code: 222804

Brand: Deltalab

Deltalab's Tamper Evident Bucket (5000ml) is a disposable container made of autoclavable polypropylene. It features a tight fitting and leakproof lid.

It also has a sturdy and convenient white plastic handle for easier carrying.

It is supplied with the lid apart and made with materials suitable for alimentary use.

Specifications:

Capacity (ml)	5000
External Lid Diameter (mm)	225
Internal Mouth Diameter (mm)	211
Height Without Lid (mm)	194
Body Weight (g)	152
Material	Polypropylene
Sterility	Non-Sterile
Quantity	1

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com