

Deltalab White Screw Cap Containers (200ml)

Product Code: 309903

Brand: Deltalab

Manufactured in polystyrene with a white screw cap made of polyethylene, Deltalab's White Screw Cap Containers (200ml) are ideal for sample collection, storage, and transport. They also provide great chemical and thermal resistance.



Specifications:

Volume (ml)	200
Max. Capacity (ml)	200
Length (mm)	58
Height (mm)	95
Cap Length (mm)	62
Cap Height (mm)	22
Material(s)	Polystyrene (container) and polyethylene (cap)
Sterility	Non-Sterile
Quantity	300

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com