

Deltalab Yellow Wide Mouth Wash Bottles (500ml)

Product Code: 191638.06

Brand: Deltalab

Deltalab's Yellow Wide Mouth Wash Bottles (500ml) are translucent low density polyethylene bottles with a coloured screw cap and delivery tube. Suitable for a variety of liquid storage and sampling applications, the tops of these bottles provide enable easy identification. The wide mouth is used to provide fast and easy filling.

Features:

- Feature a one piece cap and delivery tube assembly, thus eliminating air or liquid absorption through the cap.
- Includes a wide neck for safe and easy filling.
- Ribbed caps for a better handling while opening and closing.
- Graduated in 100ml increments.

Specifications:

Capacity (ml)	500
Cap Colour	Yellow
Internal Neck (Ø mm)	38
Base (Ø mm)	75
Height (mm)	165
Material	Polyethylene
Sterility	Non-Sterile
Quantity	1



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com