

## Deltalab Yellow Wide Mouth Wash Bottles (500ml)

**Product Code:** I9I638.06

**Brand:** Deltalab

Deltalab's Yellow Wide Mouth Wash Bottles (500ml) are translucent low density polyethylene bottles with a coloured screw cap and delivery tube. Suitable for a variety of liquid storage and sampling applications, the tops of these bottles provide enable easy identification. The wide mouth is used to provide fast and easy filling.

### Features:

- Feature a one piece cap and delivery tube assembly, thus eliminating air or liquid absorption through the cap.
- Includes a wide neck for safe and easy filling.
- Ribbed caps for a better handling while opening and closing.
- Graduated in 100ml increments.

### Specifications:

<b>Capacity (ml)</b>	500
<b>Cap Colour</b>	Yellow
<b>Internal Neck (Ø mm)</b>	38
<b>Base (Ø mm)</b>	75
<b>Height (mm)</b>	165
<b>Material</b>	Polyethylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	10



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)