

Medline 5000ml Digital Mantle

Product Code: DM-607

Brand: Medline Scientific

The Medline 5000ml Digital Mantle features an integral digital temperature controller linked to a sensor probe, which is inserted directly into the flask. The temperature is shown clearly on an LED display and the desired temperature is controlled by a simple push button potentiometer.

A retort rod holder is integrated into the base of the unit to accommodate our range of accessories.



Specifications

Volume (ml)	5000
Case Construction	Steel, with chemical resistant coating
Thermal Insulation	Moulded silica powder
Heating Element	Nichrome wire
Mantle	Woven glass fibre mat
Element Temperature	Up to 400 °C (max.)
Sensor Probe	CA (K) Type, Ø 4.8mm * 200mm (1.5m lead)
Power Supply	220V AC
Temperature Range (options)	Low (up to 100°C)/High (up to 200°C)
Watts	700

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com