

## Medline 10,000ml Digital Mantle

**Product Code:** DM-608

**Brand:** Medline Scientific

The Medline 10,000ml Digital Mantle features an integral digital temperature controller linked to a sensor probe, which is inserted directly into the flask. The temperature is shown clearly on an LED display and the desired temperature is controlled by a simple push button potentiometer.

A retort rod holder is integrated into the base of the unit to accommodate our range of accessories.



## Specifications

<b>Volume (ml)</b>	10,000
<b>Case Construction</b>	Steel, with chemical resistant coating
<b>Thermal Insulation</b>	Moulded silica powder
<b>Heating Element</b>	Nichrome wire
<b>Mantle</b>	Woven glass fibre mat
<b>Element Temperature</b>	Up to 400 °C (max.)
<b>Sensor Probe</b>	CA (K) Type, Ø 4.8mm * 200mm (1.5m lead)
<b>Power Supply</b>	220V AC
<b>Temperature Range (options)</b>	Low (up to 100°C)/High (up to 200°C)
<b>Watts</b>	1200

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)