

Medline 5000ml Controlled Heating Mantle

Product Code: E-I07N

Brand: Medline Scientific

The Medline 5000ml Controlled Heating Mantle is an energy-regulated, all-purpose heating mantle.

Its time proportional double pole system allows stable control of the heating element, reducing the risk of overheating liquids. The dial controller allows accurate temperature settings. A retort rod holder is integrated into the base of the unit to accommodate our range of accessories.



Specifications

Volume (ml)	5000
Case Construction	Steel, with chemical resistant coating
Thermal Insulation	Moulded silica powder
Heating Element	Nichrome wire
Mantle	Woven glass fibre mat
Element Temperature	Up to 450 °C
Power Supply	220V AC
Watts	700

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com