

## Medline 10,000ml Controlled Heating Mantle

**Product Code:** E-108N

**Brand:** Medline Scientific

The Medline 10,000ml Controlled Heating Mantle is a an energy-regulated, all-purpose heating mantle.

Its time proportional double pole system allows stable control of the heating element, reducing the risk of overheating liquids. The dial controller allows accurate temperature settings. A retort rod holder is integrated into the base of the unit to accommodate our range of accessories.



## Specifications

<b>Volume (ml)</b>	10,000
<b>Case Construction</b>	Steel, with chemical resistant coating
<b>Thermal Insulation</b>	Moulded silica powder
<b>Heating Element</b>	Nichrome wire
<b>Mantle</b>	Woven glass fibre mat
<b>Element Temperature</b>	Up to 450 °C
<b>Power Supply</b>	220V AC
<b>Watts</b>	1200

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)