

## Medline 5000ml Heating Stirrer Mantle

**Product Code:** ES-307N

**Brand:** Medline Scientific

The Medline 5000ml Heating Stirrer Mantle features an integral magnetic stirrer, which is independently controlled.

The stirrer speed may be varied from 100 to 1000rpm and the temperature is regulated by an infinity bimetallic thermostat. A retort rod holder is integrated into the base of the unit to accommodate our range of accessories.



## Specifications

<b>Volume (ml)</b>	5000
<b>Case Construction</b>	Steel, with chemical resistant coating
<b>Thermal Insulation</b>	Moulded silica powder
<b>Heating Element</b>	Nichrome wire
<b>Mantle</b>	Woven glass fibre mat
<b>Element Temperature</b>	Up to 450 °C
<b>Speed Range</b>	10 - 1000rpm
<b>Power Supply</b>	220V AC
<b>Watts</b>	700

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)