

Medline 20,000ml Heating Stirrer Mantle

Product Code: ES-309N

Brand: Medline Scientific

The Medline 20,000ml Heating Stirrer Mantle features an integral magnetic stirrer, which is independently controlled.

The stirrer speed may be varied from 100 to 1000rpm and the temperature is regulated by an infinity bimetallic thermostat. A retort rod holder is integrated into the base of the unit to accommodate our range of accessories.



Specifications

Volume (ml)	20,000
Case Construction	Steel, with chemical resistant coating
Thermal Insulation	Moulded silica powder
Heating Element	Nichrome wire
Mantle	Woven glass fibre mat
Element Temperature	Up to 450 °C
Speed Range	10 - 1000rpm
Power Supply	220V AC
Watts	1800

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com