


Deltalab Stainless Steel Rack (for 50 x 12mm Tubes)

Product Code: F-210

Brand: Deltalab

Deltalab's Stainless Steel Rack (for 50 x 12mm Tubes) can hold 50 glass or plastic test tubes with a diameter of 12mm. Each laboratory rack has holes with a diameter of 13mm.



Specifications:

Dimensions (l x w x h)	155 x 80 x 35mm
Containers Capacity	50 (10 x 5)
Material	Stainless Steel
Sterility	Non-Sterile
Quantity	1

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com