

Deltalab Stainless Steel Rack for 50 Tubes (12mm)

Product Code: F-211

Brand: Deltalab

Deltalab's Stainless Steel Rack for 50 Tubes (12mm) can hold 50 glass or plastic tubes with a diameter of 12mm. This laboratory rack has holes with a diameter of 13mm.

Specifications:

Dimensions (l x w x h)	150 x 80 x 70mm
Tubes Capacity	50 (10 x 5)
Material	Stainless Steel
Quantity	1

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com