


## Deltalab Stainless Steel Rack for 50 Tubes (12mm)

**Product Code:** F-211

**Brand:** Deltalab

Deltalab's Stainless Steel Rack for 50 Tubes (12mm) can hold 50 glass or plastic tubes with a diameter of 12mm. This laboratory rack has holes with a diameter of 13mm.



### Specifications:

<b>Dimensions (l x w x h)</b>	150 x 80 x 70mm
<b>Tubes Capacity</b>	50 (10 x 5)
<b>Material</b>	Stainless Steel
<b>Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)