

Deltalab Styrofoam Microtube Rack with Lid (for 50 Microtubes)

Product Code: G-050

Brand: Deltalab

Ideal for long-term storage of microtubes in laboratory freezers, Deltalab's Styrofoam Microtube Rack with Lid (for 50 Microtubes) features alphanumeric graduations to enable accurate identification of the tubes.

Specifications:

Description	50 wells rack
Suitable Microtube Sizes (ml)	1.5 - 2.0
Dimensions (mm)	210 x 107 x 65
Material	Expandable Polystyrene
Sterility	Non-Sterile
Quantity	1

