

## Deltalab Styrofoam Microtube Rack with Lid (for 100 Microtubes)

**Product Code:** G-100

**Brand:** Deltalab

Ideal for long-term storage of microtubes in laboratory freezers, Deltalab's Styrofoam Microtube Rack with Lid (for 100 Microtubes) features alphanumeric graduations to enable accurate identification of the tubes.

### Specifications:

<b>Description</b>	100 wells rack
<b>Suitable Microtube Sizes (ml)</b>	1.5 - 2.0
<b>Dimensions (mm)</b>	335 x 95 x 68
<b>Material</b>	Expandable Polystyrene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1

