

A&D 152g x 0.1mg Analytical Balance (HR-I50A)

Product Code: HR-I50A

Brand: AND Instruments

The A&D 152g x 0.1mg Analytical Balance (HR-I50A) is an entry-level compact analytical balance with a shatterproof, anti-static draft shield. The rotary doors also give this balance the small footprint ever in an A&D analytical balance. It features fast response times of just 2 seconds, auto on/off buttons, RS-232 interface (USB/LAN are optional) and a useful shock absorber function.

Features

- Robust construction, strong protection against dust and water
- Shock absorber function
- Auto power ON/OFF button
- B5 size footprint for use in narrow space
- Reverse backlit LCD display for easy to read weighing results
- SCF - Statistical calculation function
- Standard RS-232C
- Fast compact SHS sensor, response time of about 2 seconds
- Quick USB interface - no driver installation required (optional)
- LAN - Ethernet interface with WinCT-Plus software (optional)
- Easily removable large breeze break
- Built-in rechargeable NiMH battery pack (optional)



Specifications

Capacity	152g
Resolution	0.1mg
Pan Size	Ø 90mm

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com