

## A&D 252g x 0.1mg Analytical Balance (HR-250AZ)

**Product Code:** HR-250AZ

**Brand:** A&D Instruments

The A&D 252g x 0.1mg Analytical Balance (HR-250AZ) is an entry-level compact analytical balance with a shatterproof, anti-static draft shield. The rotary doors also give this balance the small footprint ever in an A&D analytical balance. Behind its simple design, the balance utilises proven A&D technology that delivers superior performance to all users who wish to do superior work.

### Features:

- Robust construction, strong protection against dust and water.
- Built-in internal calibration at the touch of a button.
- Shock absorber function.
- Auto power ON/OFF button.
- B5 size footprint for use in narrow space.
- Reverse backlit LCD display for easy to read weighing results.
- SCF - Statistical calculation function
- Standard RS-232C.
- Fast compact SHS sensor, response time of about 2 seconds.
- Quick USB interface - no driver installation required (optional).
- LAN - Ethernet interface with WinCT-Plus software (optional).
- Easily removable large breeze break.
- Built-in rechargeable NiMH battery pack (optional).



### Specifications:

<b>Capacity</b>	252g
<b>Resolution</b>	0.1mg
<b>Pan Size</b>	Ø 90mm

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)