

A&D 62g/252g x 0.1mg/1mg Analytical Balance (HR-25IAZ)

Product Code: HR-25IAZ

Brand: A&D Instruments

The A&D 62/252g x 0.1mg/1mg Analytical Balance (HR-25IAZ) is an entry-level compact analytical balance with a shatterproof, anti-static draft shield. The rotary doors also give this balance the small footprint ever in an A&D analytical balance. Behind its simple design, the balance utilises proven A&D technology that delivers superior performance to all users who wish to do superior work.

Features:

- Robust construction, strong protection against dust and water.
- Built-in internal calibration at the touch of a button.
- Shock absorber function.
- Auto power ON/OFF button.
- B5 size footprint for use in narrow space.
- Reverse backlit LCD display for easy to read weighing results.
- SCF - Statistical calculation function
- Standard RS-232C.
- Fast compact SHS sensor, response time of about 2 seconds.
- Quick USB interface - no driver installation required (optional).
- LAN - Ethernet interface with WinCT-Plus software (optional).
- Easily removable large breeze break.
- Built-in rechargeable NiMH battery pack (optional).

Specifications:

Capacity	62g/252g
Resolution	0.1mg/1mg
Pan Size	Ø 90mm



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com