


Deltalab Stainless Steel Rack for 24 Containers (50mm)

Product Code: I-240

Brand: Deltalab

Deltalab's Stainless Steel Rack for 24 Containers (50mm) can hold up to 24 containers and has 50mm square holes, along with a strong woven mesh base.



Specifications:

Dimensions (l x w x h)	315 x 210 x 45mm
Containers Capacity	24 (6 x 4)
Material	Stainless Steel
Sterility	Non-Sterile
Quantity	1

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com