

Jeio Tech 205 Litre Economy Incubator (IB-21E)

Product Code: 26035

Brand: Jeio Tech

Ideal for use in various scientific fields (e.g. life sciences, pharmaceuticals, clinical science, chemistry, biology), the Jeio Tech 205 Litre Economy Incubator (IB-21E) is an affordable yet exceptionally reliable incubator featuring film type heaters attached to the inner chamber's external walls, making the incubator free from noise, vibration and mechanical troubles.

Features:

Performance:

- Temperature range from ambient +3°C to 60°C.
- Microprocessor PID control/temperature calibration/automatic tuning.
- Vibration-free design ideal for yeast cultivation. Free from noise, vibration, and mechanical troubles by film type heating method.

Convenience:

- Intuitive control panel with bright LED display (0.1°C resolution).
- Dual wait on/off timer modes (1 min. to 99 hr. 59 min.).
- Convenient preset function for three most commonly used temperature settings.
- Corrosion-resistant incoloy sheath heating element and stainless steel interior.
- Rounded inner chamber corners for easy cleaning and optimal air-flow.
- Stackable up to 2 units using the optional stacking kit.

Safety:

- Automatic run after power interruption.
- Keypad lock function prevents accidental changes during operation.
- Over-temperature and over-current protection.
- Audible and visible open door alarms.



Specifications:

Chamber Volume (litres)	205	
Temperature Range (°C)	Amb. +3 to 60	
Temperature Fluctuation at 37°C	± 0.1	
Temperature Variation at 37°C	± 1.1	
Dimensions (w x d x h)	Internal	632 x 515 x 630mm
	External	740 x 630 x 910mm
No. of Wire Shelves (std./max.)	2/14	
Net Weight (kg)	57	
Electrical Requirements (230V, 50/60Hz)	2.4A	

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com