

Lec 85 Litre Under Counter Laboratory Freezer (ISU111)

Product Code: LECISU111

Brand: Lec Medical

The Lec 85 Litre Under Counter Laboratory Freezer (ISU111) is an under bench freezer that can also be placed on top of the bench. It comes with four drawers and defrost is manual. If there is no under bench space available, this lab freezer can be situated on the bench top.



Features

- **Spark-Free Interior:** Offers safe and reliable storage of volatile liquids, reagents, chemicals, poisons, and samples with a spark-free interior that reduces the risk of internal explosions
- **Lock:** Security is essential in laboratories, so this freezer comes complete with a door lock and two keys, ensuring that only authorised individuals are able to access it
- Adjustable thermostat
- Four drawers

Specifications

Net Capacity (L)	85
Temperature Range (°C)	-18 to -25
Temperature Display	No
Door Lock	Yes
External Dimensions (W x D x H)	545 x 554 x 845mm
Shelves/Drawers	0/4
Defrost Capability	Manual defrost
Spark-Free	Yes

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com