

## Jeio Tech 100ml Flask Support (FT 010 13)

**Product Code:** AAA23706

**Brand:** Jeio Tech

The Jeio Tech 100ml Flask Support (FT 010 13) has a polypropylene body that is chemically resistant to various organic solvents.

With a convenient grip for easy carrying, this multi-volume flask holder is concave recessed and secure/safe. They also have a noticeable name tag for classification.



### Specifications:

|                               |                   |
|-------------------------------|-------------------|
| <b>Heat Resistance</b>        | -30°C to 121°C    |
| <b>Capacity (ml)</b>          | 100               |
| <b>Row x Column</b>           | 1 x 3             |
| <b>Hole (d x h)</b>           | 70 x 91mm         |
| <b>Dimensions (w x d x h)</b> | 310 x 100 x 100mm |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)