

Jeio Tech 100ml Flask Support (FT 010 34)

Product Code: AAA23701

Brand: Jeio Tech

The Jeio Tech 100ml Flask Support (FT 010 34) has a polypropylene body that is chemically resistant to various organic solvents.

With a convenient grip for easy carrying, this multi-volume flask holder is concave recessed and secure/safe. They also have a noticeable name tag for classification.



Specifications:

Heat Resistance	-30°C to 121°C
Capacity (ml)	100
Row x Column	3 x 4
Hole (d x h)	70 x 91mm
Dimensions (w x d x h)	384 x 265 x 98mm

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com