

Jeio Tech 4.8 Litre Cold Trap Bath (CTB-10)

Product Code: AAH62012K

Brand: Jeio Tech

Jeio Tech's 4.8 Litre Cold Trap Bath (CTB-10) provides effective capture of moisture or condensate with freezing temperatures down to -40°C . It is also useful as a quick cooling bath.

With its built-in dissolving system, ice or solid matter stuck to the bath can be easily removed.

Features:

Performance:

- Temperature control down to ambient -40°C .
- Built-in dissolving system to easily remove ice or solid matter stuck on the bath.

Convenience:

- LED temperature display.
- Easy installation and maintenance of glass traps (two glass traps provided).
- Quick-release front grill for easy cleaning of the condenser.
- Compact size and silent operation.
- Durable stainless steel bath.
- Transparent acrylic lides on top of the glass trap to check condensate level in the trap.

Safety:

- Complete safety protection system with warning alarms. Over-temperature, over-current protection.



Specifications:

Temperature Range (°C)		Down to -40
Material(s)	Bath (External)	Stainless steel, 1.2t, double painted and baked
	Bath (Up Desk)	Stainless steel, 1.5t
	Insulation	Polyurethane foam (30mm)
	Refrigerator	1/2 HP, R-404A (LBP)
Dimensions	Volume (litres)	4.8
	Bath (Ø x h)	210 x 204
	Glass Trap (Ø x h)	60 x 262
	External (w x d x h)	360 x 465 x 555mm
	Net Weight (kg)	46
Electrical Requirements (230V, 50Hz)		4.5A

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com