

## Jeio Tech 9.5 Litre Electric Aspirator (VE-II)

**Product Code:** AAH71015K

**Brand:** Jeio Tech

With its built-in circulating pump and water tank, the clean and efficient Jeio Tech 9.5 Litre Electric Aspirator (VE-II) is a portable dual-channel water-jet aspirator pump suitable for rotary vacuum evaporation, decompression distillation, vacuum drying, and vacuum filtering.

Its built-in check valve prevents a backflow of water entering the aspirator pump. The water tank features a drain port and a spigot in order for the water to be easily changed.

Constructed with corrosion-resistant wetted parts with 304 stainless steel, polypropylene, silicone, or nickel-coated brass.



### Specifications:

<b>Bath Volume (litres)</b>	9.5
<b>Max. Vacuum (Mpa/ mmHg)</b>	0.0906/680
<b>Free-Air Capacity</b>	18L/minute (0.6cfm) x 2ea
<b>Material (bath/aspirator)</b>	Polypropylene/Nickel coated brass x 2ea
<b>Bath Interior (w x l, d)</b>	298 x 310, 227mm
<b>Exterior (w x l x h)</b>	330 x 265 x 390mm
<b>Inhale Nozzle</b>	External Ø 9.5mm (0.4") x 2ea
<b>Drain Cook, Over-Flow Nozzle (Ø mm)</b>	17mm
<b>Net Weight (kg)</b>	5.7
<b>Electrical Requirements</b>	230V, 50Hz/0.7A

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)