

Jeio Tech Open-Ring Cover with One Opening for BW-05 Water Baths (110mm)

Product Code: AAA45524

Brand: Jeio Tech


The Jeio Tech Open-Ring Cover with One Opening for BW-05 Water Baths (110mm) allows beakers or erlenmeyer flasks to be placed directly onto the perforated stainless steel half-shelf adjuster or water tank boom.

By adjusting the ring size for each glassware, liquid loss due to evaporation is minimised.

Suitable for Jeio Tech 3.5 Litre Water Bath (BW-05H), Jeio Tech 3.5 Litre Water Bath (BW-05B), and Jeio Tech 3.5 Litre Water Bath (BW3-05G).

Specifications:

No. of Openings	1
Diameter (mm)	110
Suitable For	BW-05H/BW-05B/BW3-05G



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com