


Jeio Tech Perforated Shelf for ILP-12 (305 x 299mm)

Product Code: AAA24502

Brand: Jeio Tech

The Jeio Tech Perforated Shelf for ILP-12 (305 x 299mm) is made of stainless steel and can be readily removed without using tools for easy cleaning.

Suitable for the Jeio Tech 48 Litre Compact Low Temperature Incubator (ILP-12).



Specifications:

Material	Stainless Steel
Dimensions (w x l)	305 x 299mm
Suitable For	ILP-12
Max. Shelves Per Incubator	10

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com