

Jeio Tech Stainless Steel Flask Clamp (100ml)

Product Code: 23551

Brand: Jeio Tech

The Jeio Tech Stainless Steel Flask Clamp (100ml) is designed for 100ml Erlenmeyer flasks when using laboratory incubators and shakers.



Specifications:

Erlenmeyer Flask Size (ml)		100
Max. Mountable Flask Clamps	IST-3075(R)	24
	IST-4075(R)	39
	ISF-7100(R)	68
	ISF-7200(R)	
	ISS-7100(R)	
	ISS-7200(R)	
	ISS-3075(R)	24
	ISS-4075(R)	39
	OS-2000	16
	OS-3000	28

Unit 3, Tower Business Park
 Warpsgrove Lane
 Chalgrove
 Oxfordshire
 OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com